

## Create cell suspension from Human Tonsils

### I – Principle

---

This SOP describes the process of isolating cells from tonsils, followed by storage in liquid nitrogen.

### II – Safety Overview

---

- It is the responsibility of all personnel handling human blood samples to ensure they are vaccinated against Hepatitis B and are trained to treat all biological samples as potentially infectious.
- The Safety and Ethics Committee must approve all work before commencement.
- Standard personal protective equipment (PPE) must be worn when performing laboratory work unless otherwise stated. Standard PPE includes.
  - o Long-sleeved lab gown, enclosed footwear, safety glasses and gloves.
- Please be aware of the risks involved in working with human samples.
  - o Care must be taken when handling blood and purified cell populations obtained from clinical samples.
- Class II Biosafety Cabinet must be used when handling unfixed blood.
  - o Cabinet must be sealed and exposed to UV light for 20 min prior to use.
  - o Cabinet surface must then be wiped down with 80% (v/v) ethanol.

Any items for use in the cabinet (including pipettes) must be decontaminated with 80% (v/v) ethanol.

### III – Equipment, Reagents and Consumables

---

#### Equipment

- Finnpiquette™ F2: P2, P20, P200, and P1000
- Finnpiquette™ F2: Multichannel pipette 30-300uL
- Eppendorf Centrifuge 5910
- Class 2 Biosafety Cabinet
- CO2 incubator (at 37C)
- Mr. Frosty containers
- CellDrop – Automatic cell Counter
- Water Bath (at 37C)
- Liquid Nitrogen tank
- Esky – for transport of samples between facilities
- Fridge 4C
- Freezer -20C
- Freezer -80C

#### Reagents/Chemicals

- Sodium Pyruvate (Thermo Fisher -Gibco # 11360070)

- DeNovix Acridine Orange / Propidium Iodide (Cat# CD-AO-PI1.5)
- DMSO (Merck)
- FBS (Thermo Fisher – Gibco # 10099141)
- Ham's F12 media (Thermo Fisher – Cat# 11765054)
- Insulin-Transferrin-Selenium (Thermo Fisher – Gibco # 41400045)
- NEAA (Non-essential Amino Acids) (Thermo Fisher - Gibco # 11140050)
- Normocin (Jomar – INV-ant-nr-2)
- PBS- 1X Filtered Phosphate Buffer Saline
- Pen/Strep (Thermo Fisher -15140122)
- RPMI 1640 with L-glutamine (Thermo Fisher - Gibco # 11875119)

## Consumables

- Tips (various manufacturers) 10, 20, 200 and 1000 µL for use with Finnpiquette™ F2 pipettes (listed above)
- 1.5 ml polypropylene tubes (Various manufacturers)
- 50 ml polypropylene tubes (Various manufacturers)
- Petri Dishes (Various manufacturers) – 10 ml
- 1ml disposable transfer pipette (Various manufacturers)
- Cryovials (Corning)
- Scissors/tweezers/scalp blades
- Syringe Terumo 10 ml

## Notes:

### Antimicrobial bath (AB) contains:

- Ham's F12 media. Contains L-glutamine and phenol red
- Normocin: 50 mg/mL, use at 100 ug/mL (1/500)
- Pen/Strep Contains 10,000 units/mL of penicillin & 10,000 µg/mL of streptomycin. Use at 50 U/uL (1/200).

Use within 1 month.

### Tonsil Processing Media (TPM) contains:

- RMPI-1640. Contains L-glutamine, phenol red, and 2.0g/L sodium bicarb. Use 450 ml.
- Heat inactivated FBS, use at 10%.
- Normocin 50 mg/mL at 100 ug/mL (1/500).
- Pen/Strep Contains 10,000 units/mL of penicillin & 10,000 µg/mL of streptomycin. Use at 50 U/uL ug/ml (1/200).
- 100 x NEAA (Non-essential Amino Acids), use at 1X.
- 100 x Insulin-Transferrin-Selenium, use at 1X.
- 100 mM Sodium Pyruvate, use at 1mM.

This media is for processing of the tonsils; it is NOT SUITABLE once growing organoids. This is because the RPMI contains L-glutamine, which degrades into ammonia overtime.

## IV – SOPS

## Preparation

Process	Steps
Prepare Reagents	<ol style="list-style-type: none"> <li>1. Thaw any reagents and place at 4°C</li> <li>2. Cool centrifuge to 4°C</li> <li>3. Place Mr Frosty containers and PBS into the fridge</li> <li>4. Prepare antimicrobial bath (AB) and tonsil processing media (TPM), and keep at 4°C.</li> </ol>
Antimicrobial Bath	<ol style="list-style-type: none"> <li>5. <u>NOTE:</u> All samples should be placed into an antimicrobial bath upon receipt</li> <li>6. Aspirate PBS from 50 mL falcon tube</li> <li>7. Add 10 mL AB to tube.</li> <li>8. Incubate for 1 hour at 4 °C.</li> <li>9. During incubation, print labels for cryovials</li> </ol>

## Dissecting into single cell suspensions

Process	Steps
Tonsil's preparation	<p><u>NOTE:</u> From here, process 2 tonsils at a time. Leave other samples in AB bath at 4°C</p> <ol style="list-style-type: none"> <li>1. Aspirate AB from 50 mL tube, replace with 15 mL PBS to rinse.</li> <li>2. Place a 10 mL petri dish onto ice and fill with 15 mL of TPM.</li> <li>3. Place lid face-up inside the hood</li> <li>4. Use this to rest sterile tools later.</li> <li>5. Place tonsils into petri dish to begin dissection.</li> </ol>
Tonsil's dissection	<p><u>NOTE:</u> While dissecting, keep pieces submerged in media as much as possible</p> <ol style="list-style-type: none"> <li>6. Remove any blood clots or fibroid tissue and discard.</li> <li>7. Record observations of any blood clots, fibrous or cauterised tissue</li> <li>8. Dissect the tonsil into &lt;0.5 cm pieces.</li> </ol>
Make cell suspension	<p>Two dish method:</p> <ol style="list-style-type: none"> <li>9. Place a second 10 mL dish containing 15 mL TPM on bench (inside the hood – but not on ice)</li> <li>10. Place four 100 um cell strainers into the dish, and transfer tonsil pieces into these</li> </ol>

	<ol style="list-style-type: none"> <li>11. Using a 10 mL syringe plunger, push tonsil fragments through the cell strainers to make a cell suspension.</li> <li>12. Place a new 100um cell strainer onto a 50 mL tube.</li> <li>13. Using a 10 mL stripette, transfer cell suspension from the petri dish over the 100-um filter and into the 50 mL tube</li> <li>14. Remove all cell strainers from the petri dish.</li> <li>15. Rinse the petri dish with 10 mLs of TPM, then pipette this volume over the cell strainer and into the 50 mL tube.</li> <li>16. Rinse the 100um strainer with a further 5 mL TPM, collect in the 50 mL tube.</li> <li>17. Incubate the resulting cell suspension at 4°C for 10 minutes.</li> <li>18. The stroma settles to bottom of tube.</li> <li>19. Gently transfer cell suspension onto a new cell strainer and into to a new 50 mL tube, leaving behind the settled stroma</li> <li>20. Top up tube to 45 mLs with TPM</li> <li>21. Rinse cell strainer in the process</li> <li>22. Centrifuge the cell suspension at 300 g for 10 minutes, acceleration 6, brake 6, at 4°C</li> <li>23. Aspirate supernatant</li> <li>24. Resuspend in 20-30 mL TPM.</li> <li>25. Prepare a 1/100 dilution of the sample.</li> <li>26. 10 uL samples + 990 uL TPM</li> <li>27. Count 1/100 diluted sample on cell drop.</li> <li>28. Add 10 uL of 1/100 dilution to 10 uL AO/PI</li> <li>29. Use primary cell AO/PI program with Boyle lab settings.</li> <li>30. Record on L: drive Count Record</li> </ol>
--	---

### Storing cell suspensions

Process	Steps
Sample preparation for long-term storage	<ol style="list-style-type: none"> <li>1. Centrifuge the cell suspension at 300 g for 10 minutes, acceleration 6, brake 6, at 4°C.</li> <li>2. During the spin, prepare 10% DMSO in FBS</li> <li>3. Aspirate supernatant, then resuspend cells to 100 e6 cells/mL in 10% DMSO/FBS</li> <li>4. Aliquot into "SPL" brand cryovials, 1 mL/vial</li> <li>5. Place cryovials into Mr Frosty, then place Mr Frosty into -80 °C freezer overnight</li> <li>6. Transfer samples to vapour phase cryogenic storage the next morning.</li> </ol>

